

## Frizzled 9 / CD349 Rabbit mAb

Catalog No: #56432

Package Size: #56432-1 50ul #56432-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Frizzled 9 / CD349 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	Frizzled 9 / CD349 Antibody detects endogenous levels of total Frizzled 9 / CD349
Immunogen Description	A synthesized peptide derived from human Frizzled 9 / CD349
Other Names	CD349; frizzled-9; Fz-9; FZD3; FzE6
Accession No.	Uniprot:O00144
Calculated MW	64kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB: 1:500~1:2000

## Images



Western blot analysis of Frizzled 9 expression in Human Seminoma lysate.

## Product Description

Receptor for Wnt proteins. Most of frizzled receptors are coupled to the beta-catenin canonical signaling pathway, which leads to the activation of disheveled proteins, inhibition of GSK-3 kinase, nuclear accumulation of beta-catenin and activation of Wnt target genes.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.