

## WIPI1 Rabbit mAb

Catalog No: #56469

Package Size: #56469-1 50ul #56469-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	WIPI1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse Rat
Specificity	WIPI1 Antibody detects endogenous levels of total WIPI1
Immunogen Description	A synthesized peptide derived from human WIPI1
Other Names	ATG18; ATG18A; WIPI 1; WIPI-1; WIPI49;
Accession No.	Uniprot:Q5MNZ9
Calculated MW	48kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000IHC:1:50~1:100

## Images



Western blot analysis of WIPI1 expression in A375 cell lysate.

## Product Description

May play a role in autophagy. May regulate the trafficking of proteins involved in the mannose-6-phosphate receptor (MPR) recycling pathway.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.