

HNF1 beta Rabbit mAb

Catalog No: #56470



Package Size: #56470-1 50ul #56470-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

| | |
|-----------------------|--|
| Product Name | HNF1 beta Rabbit mAb |
| Host Species | Rabbit |
| Clonality | Monoclonal |
| Isotype | Rabbit IgG |
| Purification | Affinity-chromatography |
| Applications | WB |
| Species Reactivity | Human |
| Specificity | HNF1 beta Antibody detects endogenous levels of total HNF1 beta |
| Immunogen Description | A synthesized peptide derived from human HNF1 beta |
| Other Names | FJHN; HNF 1B; HNF 2; HNF1 beta; HNF1B; HNF2; HPC11; LFB3; MODY5; TCF 2; TCF2; vHNF1; |
| Accession No. | Uniprot:P35680 |
| Calculated MW | 70kDa |
| Formulation | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage | Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle. |

Application Details

WB:1:500~1:2000

Images



Western blot analysis of HNF1 beta expression in SK-OV-3 cell lysate.

Product Description

Transcription factor, probably binds to the inverted palindrome 5'-GTTAATNATTAAC-3'.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.