

Brachyury Rabbit mAb

Catalog No: #56481



Package Size: #56481-1 50ul #56481-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Brachyury Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Mouse, Rat, Human
Specificity	Brachyury Antibody detects endogenous levels of total Brachyury
Immunogen Description	A synthesized peptide derived from human Brachyury
Other Names	Brachyury homolog; Bry; Protein T; SAVA; TFT; T-box transcription factor T; TBXT; Brachyury protein;
Accession No.	Uniprot:O15178
Calculated MW	47,41kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:100~1:500ICC/IF:1:100~1:500IP:1:100FC:1:100

Images



Western blot analysis of Brachyury expression in MUG-Chor1 (human sacral bone chordoma) cell lysate.

Product Description

Involved in the transcriptional regulation of genes required for mesoderm formation and differentiation. Binds to a palindromic site (called T site) and activates gene transcription when bound to such a site.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.