

## Vav1/2/3 Rabbit mAb

Catalog No: #56485

Package Size: #56485-1 50ul #56485-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Vav1/2/3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP FC
Species Reactivity	Human
Specificity	Vav1/2/3 Antibody detects endogenous levels of total Vav1/2/3
Immunogen Description	A synthesized peptide derived from human Vav1/2/3
Other Names	p95; p95vav; Protein vav 1; Protein vav 2; Protein vav 3; VAV; vav-T;
Accession No.	Uniprot:P15498/P52735/Q9UKW4
Calculated MW	98kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000 ICC/IF:1:50~1:200 IP:1:50 FC:1:50

## Images



Western blot analysis of Vav1/2/3 expression in Jurkat cell lysate.

## Product Description

The Vav family are Rho/Rac guanosine nucleotide exchange factors (GEFs), consisting of three members in mammalian cells (Vav, Vav2, Vav3) and one in nematodes (CelVav). First discovered based on its transforming properties, Vav is expressed mainly in hematopoietic cells and a few non-hematopoietic tissues, such as the pancreas and tooth enamels.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.