

## RSK4 Rabbit mAb

Catalog No: #56550

Package Size: #56550-1 50ul #56550-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	RSK4 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human
Specificity	RSK4 Antibody detects endogenous levels of total RSK4
Immunogen Description	A synthesized peptide derived from human RSK4
Other Names	p90RSK6; pp90RSK4; Ribosomal S6 kinase 4; rps6ka6; RSK 4;
Accession No.	Uniprot:Q9UK32
Calculated MW	84kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200FC:1:50

## Images



Western blot analysis of RSK4 expression in SH-SY5Y cell lysate.

## Product Description

Serine/threonine kinase that may play a role in mediating the growth-factor and stress induced activation of the transcription factor CREB.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.