

## Apolipoprotein CI Rabbit mAb

Catalog No: #56566



Package Size: #56566-1 50ul #56566-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Apolipoprotein CI Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human
Specificity	Apolipoprotein CI Antibody detects endogenous levels of total Apolipoprotein CI
Immunogen Description	A synthesized peptide derived from human Apolipoprotein CI
Other Names	APO C1; Apo CI; APOC 1; ApoC I; APOC1; APOC1B; ApolipoproteinCI;
Accession No.	Uniprot:P02654
Calculated MW	9kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:1000IHC:1:50~1:200ICC/IF:1:50~1:200IP:1:50FC:1:50

## Images



Western blot analysis of Apolipoprotein CI expression in Human plasma lysate.

## Product Description

Appears to modulate the interaction of APOE with beta-migrating VLDL and inhibit binding of beta-VLDL to the LDL receptor-related protein. Binds free fatty acids and reduces their intracellular esterification.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.