

CD97 Rabbit mAb

Catalog No: #56568

Package Size: #56568-1 50ul #56568-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	CD97 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC FC
Species Reactivity	Human Rat
Specificity	CD97 Antibody detects endogenous levels of total CD97
Immunogen Description	A synthesized peptide derived from human CD97
Other Names	ADGRE5; CD97; TM7LN1;
Accession No.	Uniprot:P48960
Calculated MW	55/28kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


Application Details

WB:1:500~1:2000

IHC:1:50~1:200

FC:1:50

Images



Western blot analysis of CD97 expression in human tonsil cell lysate.

Product Description

Receptor potentially involved in both adhesion and signaling processes early after leukocyte activation. Plays an essential role in leukocyte migration.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.