

TNNI1 Rabbit mAb

Catalog No: #56589

Package Size: #56589-1 50ul #56589-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TNNI1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human
Specificity	TNNI1 Antibody detects endogenous levels of total TNNI1
Immunogen Description	A synthesized peptide derived from human TNNI1
Other Names	SSTNI; TNN1; Tnni1; Troponin I;
Accession No.	Uniprot:P19237
Calculated MW	22kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200IP:1:50

Images



Western blot analysis of TNNI1 expression in Human fetal heart lysate.

Product Description

Troponin I is the inhibitory subunit of troponin, the thin filament regulatory complex which confers calcium-sensitivity to striated muscle actomyosin ATPase activity.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.