

## Nogo Receptor Rabbit mAb

Catalog No: #56600

Package Size: #56600-1 50ul #56600-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	Nogo Receptor Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human Mouse
Specificity	Nogo Receptor Antibody detects endogenous levels of total Nogo Receptor
Immunogen Description	A synthesized peptide derived from human Nogo Receptor
Other Names	NgR; Nogo 66 receptor; Nogo receptor; Ngor; Reticulon 4 receptor; Rtn4r;
Accession No.	Uniprot:Q9BZR6
Calculated MW	66kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:1000~1:5000IP:1:50

## Images



Western blot analysis of Nogo Receptor expression in HepG2 cell lysate.

## Product Description

Nogo Receptor mediates axonal growth inhibition and may play a role in regulating axonal regeneration and plasticity in the adult central nervous system. It is also a receptor for Nogo A (RTN4), OMG and MAG.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.