

ZA2G / AZGP1 Rabbit mAb

Catalog No: #56604

Package Size: #56604-1 50ul #56604-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

| | |
|-----------------------|--|
| Product Name | ZA2G / AZGP1 Rabbit mAb |
| Host Species | Rabbit |
| Clonality | Monoclonal |
| Isotype | Rabbit IgG |
| Purification | Affinity-chromatography |
| Applications | WB IHC IP |
| Species Reactivity | Human |
| Specificity | ZA2G / AZGP1 Antibody detects endogenous levels of total ZA2G / AZGP1 |
| Immunogen Description | A synthesized peptide derived from human ZA2G / AZGP1 |
| Other Names | Azgp1; ZA2G; ZAG; ZNGP1; |
| Accession No. | Uniprot:P25311 |
| Calculated MW | 43kDa |
| Formulation | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage | Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle. |

Application Details

WB:1:500~1:2000IHC:1:50~1:200IP:1:40

Images



Western blot analysis of ZA2G / AZGP1 expression in HepG2 cell lysate.

Product Description

Stimulates lipid degradation in adipocytes and causes the extensive fat losses associated with some advanced cancers. May bind polyunsaturated fatty acids.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.