

TCF21 Rabbit mAb

Catalog No: #56612

Package Size: #56612-1 50ul #56612-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TCF21 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP FC
Species Reactivity	Human
Specificity	TCF21 Antibody detects endogenous levels of total TCF21
Immunogen Description	A synthesized peptide derived from human TCF21
Other Names	Capsulin; Epicardin; POD 1; POD1; TCF-21; tcf21;
Accession No.	Uniprot:O43680
Calculated MW	26kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IP:1:60FC:1:80

Images



Western blot analysis of TCF21 expression in 293 cell lysate.

Product Description

Involved in epithelial-mesenchymal interactions in kidney and lung morphogenesis that include epithelial differentiation and branching morphogenesis. May play a role in the specification or differentiation of one or more subsets of epicardial cell types.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.