

## Bmf Rabbit mAb

Catalog No: #56614

Package Size: #56614-1 50ul #56614-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	Bmf Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human Mouse Rat
Specificity	Bmf Antibody detects endogenous levels of total Bmf
Immunogen Description	A synthesized peptide derived from human Bmf
Other Names	Bcl2 modifying factor; Bmf;
Accession No.	Uniprot:Q96LC9
Calculated MW	24kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000FC:1:50

## Images



Western blot analysis of Bmf expression in Ramos cell lysate.

## Product Description

May play a role in apoptosis. Isoform 1 seems to be the main initiator.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.