

CDKN2C Rabbit mAb

Catalog No: #56630

Package Size: #56630-1 50ul #56630-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CDKN2C Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human Mouse Rat
Specificity	CDKN2C Antibody detects endogenous levels of total CDKN2C
Immunogen Description	A synthesized peptide derived from human CDKN2C
Other Names	Cdkn2c; CDKN6; INK4C; p18; p18 inhibits CDK4; p18 INK4c; p18 INK6; p18-INK4c; p18-INK6;
Accession No.	Uniprot:P42773
Calculated MW	15kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details


WB:1:500~1:2000

IHC:1:50~1:200

ICC/IF:1:50~1:200

IP:1:40

Images



Western blot analysis of CDKN2C expression in Ramos cell lysate.

Product Description

Interacts strongly with CDK6, weakly with CDK4. Inhibits cell growth and proliferation with a correlated dependence on endogenous retinoblastoma protein RB.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.