

MDH1 Rabbit mAb

Catalog No: #56648

Package Size: #56648-1 50ul #56648-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MDH1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP FC
Species Reactivity	Human Mouse Rat
Specificity	MDH1 Antibody detects endogenous levels of total MDH1
Immunogen Description	A synthesized peptide derived from human MDH1
Other Names	MDH s; mdh1; MDHA; MDHs; MOR2;
Accession No.	Uniprot:P40925
Calculated MW	36kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:1000~1:5000|ICC/IF:1:50~1:200|IP:1:50|FC:1:30

Images



Western blot analysis of MDH1 expression in HeLa cell lysate.

Product Description

Belongs to the LDH/MDH superfamily. MDH type 2 family.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.