

VAV3 Rabbit mAb

Catalog No: #56652

Package Size: #56652-1 50ul #56652-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	VAV3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP FC
Species Reactivity	Human Mouse
Specificity	VAV3 Antibody detects endogenous levels of total VAV3
Immunogen Description	A synthesized peptide derived from human VAV3
Other Names	RGD1565941; VAV 3; Vav3; VAV3 oncogene;
Accession No.	Uniprot:Q9UKW4
Calculated MW	98kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IC/IF:1:50~1:200IP:1:50FC:1:40.

Images



Western blot analysis of VAV3 expression in Jurkat cell lysate.

Product Description

Plays an important role in angiogenesis. Its recruitment by phosphorylated EPHA2 is critical for EFNA1-induced RAC1 GTPase activation and vascular endothelial cell migration and assembly.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.