

## Neurogenin 2 Rabbit mAb

Catalog No: #56658

Package Size: #56658-1 50ul #56658-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	Neurogenin 2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	Neurogenin 2 Antibody detects endogenous levels of total Neurogenin 2
Immunogen Description	A synthesized peptide derived from human Neurogenin 2
Other Names	Atoh4; bHLHa8; Math4A; MGC; NEUROG2; Neurogenin2; NGN2;
Accession No.	Uniprot:Q9H2A3
Calculated MW	25kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB: 1:500~1:2000

## Images



Western blot analysis of Neurogenin 2 expression in SH-SY5Y cell lysate.

## Product Description

Involved in neurogenesis.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.