

## HOXA13 Rabbit mAb

Catalog No: #56671



Package Size: #56671-1 50ul #56671-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	HOXA13 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	HOXA13 Antibody detects endogenous levels of total HOXA13
Immunogen Description	A synthesized peptide derived from human HOXA13
Other Names	Homeobox1J; HomeoboxA13; HOX 1; HOX 1J; HOX A13; HOX1; HOX1J; HOXA 13; HOXA13; Transcription factor HOXA13;
Accession No.	Uniprot:P31271
Calculated MW	42kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB: 1:500~1:2000

## Images



Western blot analysis of HOXA13 expression in 293T cell lysate.

## Product Description

Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.