

ApoD Rabbit mAb

Catalog No: #56687

Package Size: #56687-1 50ul #56687-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ApoD Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	ApoD Antibody detects endogenous levels of total ApoD
Immunogen Description	A synthesized peptide derived from human ApoD
Other Names	ApoD; ApolipoproteinD;
Accession No.	Uniprot:P05090
Calculated MW	28kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200

Images



Western blot analysis of ApoD expression in Human plasma lysate.

Product Description

APOD occurs in the macromolecular complex with lecithin-cholesterol acyltransferase. It is probably involved in the transport and binding of bilin. Appears to be able to transport a variety of ligands in a number of different contexts.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.