

ApoA4 Rabbit mAb

Catalog No: #56694

Package Size: #56694-1 50ul #56694-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ApoA4 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human
Specificity	ApoA4 Antibody detects endogenous levels of total ApoA4
Immunogen Description	A synthesized peptide derived from human ApoA4
Other Names	Apo AIV; APOA 4; ApoA IV; APOA4; Apolipoprotein A IV; Apolipoprotein A4;
Accession No.	Uniprot:P06727
Calculated MW	45kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IP:1:50

Images



Western blot analysis of ApoA4 expression in Human plasma lysate.

Product Description

APOA4 (apolipoprotein A-IV) is a component of HDL and chylomicrons. Its primary site of synthesis is the intestine, in association with lymph chylomicron particles.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.