

## UNG Rabbit mAb

Catalog No: #56696



Package Size: #56696-1 50ul #56696-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	UNG Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	UNG Antibody detects endogenous levels of total UNG
Immunogen Description	A synthesized peptide derived from human UNG
Other Names	UNG; DGU; HIGM4; HIGM5; Uracil-DNA glycosylase; UNG1; UNG2; UDG; UNG15;
Accession No.	Uniprot:P13051
Calculated MW	32kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200FC:1:50

## Images



Western blot analysis of UNG expression in HeLa cell lysate.

## Product Description

Excises uracil residues from the DNA which can arise as a result of misincorporation of dUMP residues by DNA polymerase or due to deamination of cytosine.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.