

Thyroid Peroxidase Rabbit mAb

Catalog No: #56712

Package Size: #56712-1 50ul #56712-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

| | |
|-----------------------|--|
| Product Name | Thyroid Peroxidase Rabbit mAb |
| Host Species | Rabbit |
| Clonality | Monoclonal |
| Isotype | Rabbit IgG |
| Purification | Affinity-chromatography |
| Applications | WB IHC IP |
| Species Reactivity | Human |
| Specificity | Thyroid Peroxidase Antibody detects endogenous levels of total Thyroid Peroxidase |
| Immunogen Description | A synthesized peptide derived from human Thyroid Peroxidase |
| Other Names | MSA; TDH2A; TPO; TPX; |
| Accession No. | Uniprot:P07202 |
| Calculated MW | 103kDa |
| Formulation | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage | Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle. |

Application Details

WB:1:1000~1:5000IHC:1:100~1:500IP:1:50

Images



Western blot analysis of Thyroid Peroxidase expression in Human thyroid lysate.

Product Description

Iodination and coupling of the homonogenic tyrosines in thyroglobulin to yield the thyroid hormones T(3) and T(4).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.