

IL17A Receptor Rabbit mAb

Catalog No: #56715



Package Size: #56715-1 50ul #56715-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	IL17A Receptor Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human
Specificity	IL17A Receptor Antibody detects endogenous levels of total IL17A Receptor
Immunogen Description	A synthesized peptide derived from human IL17A Receptor
Other Names	IL-17 receptor; IL-17RA; IL17RA; interleukin-17 receptor;
Accession No.	Uniprot:Q96F46
Calculated MW	160kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:1000~1:5000IP:1:50

Images



Western blot analysis of IL17A Receptor expression in Raji cell lysate.

Product Description

Receptor for IL17A and IL17F. Binds its IL17A ligand with low affinity, suggesting that additional components are involved in IL17A-induced signaling.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.