

CTRL Rabbit mAb

Catalog No: #56727

Package Size: #56727-1 50ul #56727-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CTRL Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	CTRL Antibody detects endogenous levels of total CTRL
Immunogen Description	A synthesized peptide derived from human CTRL
Other Names	CTRL; CTRL1;
Accession No.	Uniprot:P40313
Calculated MW	24,28kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000

Images



Western blot analysis of CTRL expression in Mouse pancreas lysate.

Product Description

The function of this protein remains unknown.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.