

# Thromboxane synthase Rabbit mAb

Catalog No: #56728

Package Size: #56728-1 50ul #56728-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	Thromboxane synthase Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse
Specificity	Thromboxane synthase Antibody detects endogenous levels of total Thromboxane synthase
Immunogen Description	A synthesized peptide derived from human Thromboxane synthase
Other Names	TS; TXS; CYP5; THAS; TXAS; CYP5A1; GHOSAL; BDPLT14;
Accession No.	Uniprot:P24557
Calculated MW	54kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB: 1:500~1:2000

## Images



Western blot analysis of Thromboxane synthase expression in Human fetal liver lysate.

## Product Description

This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.