

Hippocalcin Rabbit mAb

Catalog No: #56766



Package Size: #56766-1 50ul #56766-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Hippocalcin Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	Hippocalcin Antibody detects endogenous levels of total Hippocalcin
Immunogen Description	A synthesized peptide derived from human Hippocalcin
Other Names	BDR2; HPCA; P23K;
Accession No.	Uniprot:P84074
Calculated MW	26kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000

Images



Western blot analysis of Hippocalcin expression in Neuro2a cell lysate.

Product Description

May be involved in the calcium-dependent regulation of rhodopsin phosphorylation. Binds two calcium ions.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.