

PLTP Rabbit mAb

Catalog No: #56768



Package Size: #56768-1 50ul #56768-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PLTP Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Rat
Specificity	PLTP Antibody detects endogenous levels of total PLTP
Immunogen Description	A synthesized peptide derived from human PLTP
Other Names	BPIFE; HDLCQ9; OD107; Pltp;
Accession No.	Uniprot:P55058
Calculated MW	80kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000

Images



Western blot analysis of PLTP expression in HeLa cell lysate.

Product Description

Converts HDL into larger and smaller particles. May play a key role in extracellular phospholipid transport and modulation of hdl particles.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.