

PER3 Rabbit mAb

Catalog No: #56780

Package Size: #56780-1 50ul #56780-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PER3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	PER3 Antibody detects endogenous levels of total PER3
Immunogen Description	A synthesized peptide derived from human PER3
Other Names	GIG13; gPER3; hPER3; mPer3; Per3; Period 3; Period circadian clock 3; Period homolog 3; Period3;
Accession No.	Uniprot:P56645
Calculated MW	132kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000

Images



Western blot analysis of PER3 expression in K562 cell lysate.

Product Description

Component of the circadian clock mechanism which is essential for generating circadian rhythms. Can bind heme.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.