

## MFF Rabbit mAb

Catalog No: #56797

Package Size: #56797-1 50ul #56797-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	MFF Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	MFF Antibody detects endogenous levels of total MFF
Immunogen Description	A synthesized peptide derived from human MFF
Other Names	C2orf33; GL004; Mff;
Accession No.	Uniprot:Q9GZY8
Calculated MW	34kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB: 1:500~1:2000

## Images



Western blot analysis of MFF expression in Human fetal muscle lysate.

## Product Description

Plays a role in mitochondrial and peroxisomal fission.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.