

## PI 3 Kinase R4 Rabbit mAb

Catalog No: #56805

Package Size: #56805-1 50ul #56805-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	PI 3 Kinase R4 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	PI 3 Kinase R4 Antibody detects endogenous levels of total PI 3 Kinase R4
Immunogen Description	A synthesized peptide derived from human PI 3 Kinase R4
Other Names	P150; Phosphoinositide 3 kinase regulatory subunit 4; PI3 kinase p150 subunit; PIK3R4; VPS15;
Accession No.	Uniprot:Q99570
Calculated MW	150kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200FC:1:50

## Images



Western blot analysis of PI 3 Kinase R4 expression in (1) Jurkat cell lysate; (2) Molt-4 cell lysate.

## Product Description

Regulatory subunit of the PI3K complex. May regulate membrane trafficking late in the endocytic pathway.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.