

HMGN2 Rabbit mAb

Catalog No: #56813

Package Size: #56813-1 50ul #56813-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	HMGN2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	HMGN2 Antibody detects endogenous levels of total HMGN2
Immunogen Description	A synthesized peptide derived from human HMGN2
Other Names	HMG17; HMGN2;
Accession No.	Uniprot:P05204
Calculated MW	18kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200

Images



Western blot analysis of HMGN2 expression in K562 cell lysate.

Product Description

Binds to the inner side of the nucleosomal DNA thus altering the interaction between the DNA and the histone octamer. May be involved in the process which maintains transcribable genes in an unique chromatin conformation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.