

## FMO3 Rabbit mAb

Catalog No: #56822

Package Size: #56822-1 50ul #56822-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	FMO3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	FMO3 Antibody detects endogenous levels of total FMO3
Immunogen Description	A synthesized peptide derived from human FMO3
Other Names	FMO3; FMOII; TMAU;
Accession No.	Uniprot:P31513
Calculated MW	56kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200

## Images



Western blot analysis of FMO3 expression in Mouse kidney lysate.

## Product Description

Involved in the oxidative metabolism of a variety of xenobiotics such as drugs and pesticides. It N-oxygenates primary aliphatic alkylamines as well as secondary and tertiary amines.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.