

LPP Rabbit mAb

Catalog No: #56824

Package Size: #56824-1 50ul #56824-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	LPP Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	LPP Antibody detects endogenous levels of total LPP
Immunogen Description	A synthesized peptide derived from human LPP
Other Names	LIM protein; Lipoma preferred partner; lpp;
Accession No.	Uniprot:Q93052
Calculated MW	70kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200

Images



Western blot analysis of LPP expression in HeLa cell lysate.

Product Description

May play a structural role at sites of cell adhesion in maintaining cell shape and motility. In addition to these structural functions, it may also be implicated in signaling events and activation of gene transcription.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.