

## DPYD Rabbit mAb

Catalog No: #56829

Package Size: #56829-1 50ul #56829-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	DPYD Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human
Specificity	DPYD Antibody detects endogenous levels of total DPYD
Immunogen Description	A synthesized peptide derived from human DPYD
Other Names	DHP; DHPDHase; DPD; DPYD;
Accession No.	Uniprot:Q12882
Calculated MW	110kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200

## Images



Western blot analysis of DPYD expression in HeLa cell lysate.

## Product Description

Involved in pyrimidine base degradation. Catalyzes the reduction of uracil and thymine. Also involved the degradation of the chemotherapeutic drug 5-fluorouracil.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.