

HMGCS2 Rabbit mAb

Catalog No: #56833

Package Size: #56833-1 50ul #56833-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	HMGCS2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse Rat
Specificity	HMGCS2 Antibody detects endogenous levels of total HMGCS2
Immunogen Description	A synthesized peptide derived from human HMGCS2
Other Names	HMG CoA synthase; HMGCS2;
Accession No.	Uniprot:P54868
Calculated MW	50kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200

Images



Western blot analysis of HMGCS2 expression in Caco-2 cell lysate.

Product Description

This enzyme condenses acetyl-CoA with acetoacetyl-CoA to form HMG-CoA, which is the substrate for HMG-CoA reductase.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.