

MYH7B Rabbit mAb

Catalog No: #56846

Package Size: #56846-1 50ul #56846-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MYH7B Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse Rat
Specificity	MYH7B Antibody detects endogenous levels of total MYH7B
Immunogen Description	A synthesized peptide derived from human MYH7B
Other Names	Myh7b; Myosin-7B;
Accession No.	Uniprot:A7E2Y1
Calculated MW	221kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200

Images



Western blot analysis of MYH7B expression in Human fetal muscle lysate.

Product Description

Involved in muscle contraction.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.