

HOXA5 Rabbit mAb

Catalog No: #56858

Package Size: #56858-1 50ul #56858-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	HOXA5 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	HOXA5 Antibody detects endogenous levels of total HOXA5
Immunogen Description	A synthesized peptide derived from human HOXA5
Other Names	HOX1; HOX1.3; HOX1C; HOXA5;
Accession No.	Uniprot:P20719
Calculated MW	43kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000

Images



Western blot analysis of HOXA5 expression in 293T cell lysate.

Product Description

Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis. Also binds to its own promoter. Binds specifically to the motif 5'-CYYNATTA[*TG*Y]-3'.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.