

## ANPRB Rabbit mAb

Catalog No: #56859

Package Size: #56859-1 50ul #56859-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	ANPRB Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	ANPRB Antibody detects endogenous levels of total ANPRB
Immunogen Description	A synthesized peptide derived from human ANPRB
Other Names	AMDM; ANPB; ANPRB; GCB; GUC2B; GUCY2B; Npr2; NPRB; NPRBi;
Accession No.	Uniprot:P20594
Calculated MW	117kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000FC:1:50

## Images



Western blot analysis of ANPRB expression in Human fetal kidney lysate.

## Product Description

Receptor for the C-type natriuretic peptide NPPC/CNP hormone. Has guanylate cyclase activity upon binding of its ligand. May play a role in the regulation of skeletal growth.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.