

## PCBP1 Rabbit mAb

Catalog No: #56893

Package Size: #56893-1 50ul #56893-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	PCBP1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human Mouse Rat
Specificity	PCBP1 Antibody detects endogenous levels of total PCBP1
Immunogen Description	A synthesized peptide derived from human PCBP1
Other Names	Alpha-CP1; hnRNP E1; HNRPE1; HNRPX; PCBP1;
Accession No.	Uniprot:Q15365
Calculated MW	40kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000IP:1:50

## Images



Western blot analysis of PCBP1 expression in HeLa cell lysate.

## Product Description

Single-stranded nucleic acid binding protein that binds preferentially to oligo dC.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.