

Thiophosphate ester Rabbit mAb

Catalog No: #56912

Package Size: #56912-1 50ul #56912-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Thiophosphate ester Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	IP
Specificity	Thiophosphate ester Antibody detects endogenous levels of Thiophosphate ester
Immunogen Description	A synthesized peptide derived from human Thiophosphate ester
Accession No.	Uniprot:
Calculated MW	kDa
Storage	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

IP:1:1000

Product Description

Thiophosphate ester is used as fungicide.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.