

## GBP1 Rabbit mAb

Catalog No: #56930

Package Size: #56930-1 50ul #56930-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	GBP1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human
Specificity	GBP1 Antibody detects endogenous levels of total GBP1
Immunogen Description	A synthesized peptide derived from human GBP1
Other Names	GBP1; GTP binding protein 1; Guanine nucleotide binding protein 1; Guanylate binding protein 1; HuGBP1; Interferon induced guanylate binding protein 1;
Accession No.	Uniprot:P32455
Calculated MW	68kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB: 1:500~1:2000 IHC: 1:50~1:200 ICC/IF: 1:50~1:200

## Images



Western blot analysis of GBP1 expression in HeLa cell lysate.

## Product Description

Binds GTP, GDP and GMP.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.