

UPP1 Rabbit mAb

Catalog No: #56941

Package Size: #56941-1 50ul #56941-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

| | |
|-----------------------|--|
| Product Name | UPP1 Rabbit mAb |
| Host Species | Rabbit |
| Clonality | Monoclonal |
| Isotype | Rabbit IgG |
| Purification | Affinity-chromatography |
| Applications | WB |
| Species Reactivity | Human Mouse Rat |
| Specificity | UPP1 Antibody detects endogenous levels of total UPP1 |
| Immunogen Description | A synthesized peptide derived from human UPP1 |
| Other Names | UDRPASE; UPase 1; UPASE; UPP; Upp1; UrdPase 1; |
| Accession No. | Uniprot:Q16831 |
| Calculated MW | 28kDa |
| Formulation | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage | Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle. |

Application Details

WB:1:500~1:2000

Images



Western blot analysis of UPP1 expression in (1) HeLa cell lysate; (2) 293 cell lysate.

Product Description

Catalyzes the reversible phosphorylytic cleavage of uridine and deoxyuridine to uracil and ribose- or deoxyribose-1-phosphate. The produced molecules are then utilized as carbon and energy sources or in the rescue of pyrimidine bases for nucleotide synthesis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.