

ACAA2 Rabbit mAb

Catalog No: #56943

Package Size: #56943-1 50ul #56943-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ACAA2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	ACAA2 Antibody detects endogenous levels of total ACAA2
Immunogen Description	A synthesized peptide derived from human ACAA2
Other Names	Acaa2; DSAEC;
Accession No.	Uniprot:P42765
Calculated MW	45kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


Application Details

WB: 1:500~1:2000

IHC: 1:50~1:200

ICC/IF: 1:50~1:200

Images



Western blot analysis of ACAA2 expression in HeLa cell lysate.

Product Description

Abolishes BNIP3-mediated apoptosis and mitochondrial damage.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.