

UBE4A Rabbit mAb

Catalog No: #56944

Package Size: #56944-1 50ul #56944-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	UBE4A Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human
Specificity	UBE4A Antibody detects endogenous levels of total UBE4A
Immunogen Description	A synthesized peptide derived from human UBE4A
Other Names	UBE4A; UBOX2; UFD2; UFD2b;
Accession No.	Uniprot:Q14139
Calculated MW	123kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000 ICC/IF:1:50~1:200 FC:1:50

Images



Western blot analysis of UBE4A expression in 293 cell lysate.

Product Description

Binds to the ubiquitin moieties of preformed conjugates and catalyzes ubiquitin chain assembly in conjunction with E1, E2, and E3.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.