

## BRMS1 Rabbit mAb

Catalog No: #56946



Package Size: #56946-1 50ul #56946-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	BRMS1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human
Specificity	BRMS1 Antibody detects endogenous levels of total BRMS1
Immunogen Description	A synthesized peptide derived from human BRMS1
Other Names	BRMS1;
Accession No.	Uniprot:Q9HCU9
Calculated MW	28kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details


WB:1:500~1:2000

IHC:1:50~1:200

ICC/IF:1:50~1:200

IP:1:50

## Images



Western blot analysis of BRMS1 expression in HeLa cell lysate.

## Product Description

Transcriptional repressor. Down-regulates transcription activation by NF-kappa-B by promoting the deacetylation of RELA at 'Lys-310'. Promotes HDAC1 binding to promoter regions.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.