

CD58 Rabbit mAb

Catalog No: #56951

Package Size: #56951-1 50ul #56951-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD58 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP FC
Species Reactivity	Human
Specificity	CD58 Antibody detects endogenous levels of total CD58
Immunogen Description	A synthesized peptide derived from human CD58
Other Names	AG3; CD58; LFA3;
Accession No.	Uniprot:P19256
Calculated MW	50-75kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


Application Details

WB:1:500~1:2000

IP:1:50

FC:1:50

Images



Western blot analysis of CD58 expression in Raji cell lysate.

Product Description

Ligand of the T-lymphocyte CD2 glycoprotein.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.