

CD79b Rabbit mAb

Catalog No: #56964

Package Size: #56964-1 50ul #56964-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CD79b Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human
Specificity	CD79b Antibody detects endogenous levels of total CD79b
Immunogen Description	A synthesized peptide derived from human CD79b
Other Names	B29; CD79b; IGB; Igbeta;
Accession No.	Uniprot:P40259
Calculated MW	30-50kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IC/IF:1:50~1:200FC:1:50

Images



Western blot analysis of CD79b expression in Ramos cell lysate.

Product Description

Required in cooperation with CD79A for initiation of the signal transduction cascade activated by the B-cell antigen receptor complex (BCR) which leads to internalization of the complex, trafficking to late endosomes and antigen presentation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.