

## MRRF Rabbit mAb

Catalog No: #56965

Package Size: #56965-1 50ul #56965-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	MRRF Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse
Specificity	MRRF Antibody detects endogenous levels of total MRRF
Immunogen Description	A synthesized peptide derived from human MRRF
Other Names	MRFF; MTRRF; Ribosome recycling factor, mitochondrial; RRF;
Accession No.	Uniprot:Q96E11
Calculated MW	24kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000IHC:1:50~1:200

## Images



Western blot analysis of MRRF expression in HeLa cell lysate.

## Product Description

Responsible for the release of ribosomes from messenger RNA at the termination of protein biosynthesis. May increase the efficiency of translation by recycling ribosomes from one round of translation to another.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.