

## CTNNBIP1 Rabbit mAb

Catalog No: #56967

Package Size: #56967-1 50ul #56967-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	CTNNBIP1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP FC
Species Reactivity	Human
Specificity	CTNNBIP1 Antibody detects endogenous levels of total CTNNBIP1
Immunogen Description	A synthesized peptide derived from human CTNNBIP1
Other Names	CTNNBIP1; ICAT;
Accession No.	Uniprot:Q9NSA3
Calculated MW	12kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000IP:1:50FC:1:50

## Images



Western blot analysis of CTNNBIP1 expression in 293T cell lysate.

## Product Description

Prevents the interaction between CTNNB1 and TCF family members, and acts as negative regulator of the Wnt signaling pathway.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.