

HLA DMB Rabbit mAb

Catalog No: #56971

Package Size: #56971-1 50ul #56971-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	HLA DMB Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human
Specificity	HLA DMB Antibody detects endogenous levels of total HLA DMB
Immunogen Description	A synthesized peptide derived from human HLA DMB
Other Names	DMB; RING7;
Accession No.	Uniprot:P28068
Calculated MW	29kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000FC:1:50

Images



Western blot analysis of HLA DMB expression in Raji cell lysate.

Product Description

HLA-DMB belongs to the HLA class II beta chain paralogues. This class II molecule is a heterodimer consisting of an alpha (DMA) and a beta (DMB) chain, both anchored in the membrane.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.